


BRA-135F XRF Energy Dispersive General Purpose Spectrometer



BRA-135F XRF Energy Dispersive General Purpose Spectrometer



Wide range of detected chemical elements – 9F - 92U


High sensitivity due to optimal X-ray optical path

High resolution of ultra-modern SDD detector

Fundamental parameter (FP) Method for steel and alloy quantitative analysis

Built-in control PC

Patent for X-ray transparent dividing window



Completely satisfies the requirements of radiation safety.

X-ray fluorescence energy dispersive general purpose spectrometer BRA-135F allows simultaneous determination of chemical elements by characteristic energies in the 1 to 30 keV range (where elements from F till U are fitted) over a wide scope of concentrations from hundreds ppb. BRA-135F analyzes solid, powder and liquid samples, thin layer on the surface or precipitated on filters.

Operating principle


The spectrometer operating principle is based on excitation of fluorescence radiation of atoms in the substance being analyzed by radiation coming from the low-power X-ray tube. The fluorescence radiation from the sample gets into the SDD semiconductor detector where quanta of different energies are converted into electrical pulses, which amplitudes are proportional to the energy of absorbed quanta. Pulse-frequency rate with the certain amplitude is proportional to the chemical element concentration in the sample.

Method advantages

X-ray fluorescence analysis (XRF) occupies a leading position among the other methods for determination of the quantitative elemental composition of substances. XRF advantages are as follows: nondestructive measurements, multielement determination, express method, high accuracy of analysis, wide range of measured concentrations, development level of quantitative analysis theory, possibility for quantitative analysis with absence of standard samples.

High efficiency SDD detector

The silicon drift detector (SDD) with ultra-thin entrance window allows to register x-ray radiation in wide energy range on retention of energy resolution and response.



Detector scheme

Low detection limit

Owing to optimally selected materials and thickness of primary radiation filters, high transparent X-ray optical scheme, a low detection limit can be achieved for all elements to be analyzed. In the range of light elements from ^9F to ^{17}Cl low-energy radiation registration becomes possible using vacuum where the optical path of radiation passes.

Large or irregular shape samples

It is possible to measure large-size or odd-shaped samples:

- Large-size minerals and nuggets;
- Industrial articles for analysis for ROHS requirements;
- Metals and alloys incoming control;
- Analysis of liquids in the special cells or on the special filters.

Compact body and functionality

- The housing spectrometer ensure fully radiation-protection;
- Handles for carrying;
- Built-in operational computer (PC);
- LAN port for remote control and archiving of measurement results;
- LIMS integration is available;
- Easy report creation;
- Password protection and separation of access rights.

Technical data

| | |
|--|------------------------------|
| Range of detected elements | $^9\text{F} - ^{92}\text{U}$ |
| Limits of detection without preliminary enrichment, % | |
| - for elements from Na to Mg | $n \cdot 10^{-2}$ |
| - for elements from Al to Cl | 0,002 |
| - for elements from K to U | 0,0005 |
| Limit of determination at sample preconcentration (depending on chemical element), % | $1.5 \cdot 10^{-5}$ |
| Limit of determination in mid group element (liquid), g/dm ³ | $n \cdot 10^{-3}$ |
| Average time of one sample analysis, s | 100 |
| Energy resolution on MnKa line at pulses counting rate up to 10^4 s^{-1} , not more than, eV | 145 |
| Max. voltage of X-ray tube, kV | 50 |
| X-ray power, W | 10 |
| X-ray tube cooling | by air |
| Primary X-ray radiation filters, pcs | 5 |
| Number of samples installed into sample changer, up to | |
| changer #1 (Ø34 mm samples) | 15 |
| changer #2 (Ø34,36,40,44 mm samples) | 11 |
| Maximum sample size, mm | Ø 200x60 |
| Ethernet connection | Yes |
| Possibility for remote control | Yes |
| Overall dimensions (LxWxH), mm | 700x410x400 |
| Instrument weight, max, kg | 65 |
| Power | 220 V, 50 Hz |
| Power consumption, W | 500 |



Changer is able to carry solid, liquid or powder samples



Irregular shape samples

Filed of application



Methodology description

Oil analysis

For measuring purposes of trace elements Al, Ba, Ca, Cu, Fe, Mn, V, Ni, Pb, Zn, P into oil and petrochemicals appropriate methodology was developed.

Analytic complex consisting of BRA-135F and measuring methodology is capable of carrying out quantitative element analysis of petrochemicals in order to define metal trace elements and can be used to analyse exhausted motor oils of aircraft, machines, special motor vehicle in order to identify deterioration rate of engines and define applicability of technical service for it. Methodology is purchased additionally.

Besides BRA-135F could be used for testing as per ASTM D4294 and ASTM D6481. These test methods cover the measurement of sulfur, Barium, Calcium, Magnesium, Phosphorus, Zinc & Chlorine in hydrocarbons, such as lubricating used oil, diesel, naphtha, kerosene, residuals, lubricating base oils, hydraulic oils, grease, jet fuels, crude oils, gasoline (all unleaded), and other distillates. Additionally, sulfur in other products, such as M-85 and M-100, may be analyzed using this technique.

Detection limits of BRA-135F according to certified methodology (ppm):

| P | Al | Mn | Ba | Pb | V | Cu | Ni | Fe | Zn | Ca |
|-----|-----|----|----|----|---|----|----|----|----|----|
| 100 | 100 | 5 | 50 | 5 | 5 | 5 | 5 | 5 | 5 | 50 |

Cement materials analysis

For measuring purposes of mass fraction of Na, Mg, Al, Si, P, S, Cl, K, Ca, Ti, Cr, Mn, Fe, Zn, Rb, Sr into cements and cement production materials (clinkers, raw mixes) suitable methodology was developed.

Analysis methodology includes algorithms of principal components determination using X-ray energy-dispersive fluorescent spectrometers BRA-135F and is based on recommendations from GOST 5382-91, GOST R 55410-2013 (ISO 12677:2011).

There is provided remelting method into platinum crucibles (according to GOST R 55410-2013) for samples preparation in the methodology.

[illegible]

[illegible]

Marketing, Advertisement and Sales Department:

Tel.: +7 (812) 458-89-95, 458-86-48

E-mail: marketing@bv.alrosa.ru

Aftersales Service

Tel./Fax: +7 (812) 528-82-83

E-mail: quality@bv.alrosa.ru

Our partners

Bangladesh

Milestone Instruments

Telephone: +88027253468

E-mail: mstoneinst@gmail.com

India

AlfaTech Services

Telephone: +91-11-2544 6275 / 2544 6276

E-mail: alfatech@alfatechservices.com

www.alfatechservices.com

Smart Labtech PVT Limited

Telephone: +91 40-66783744, +91 40 66624394

E-mail: info@smartlabtech.net

smartlabtech.net

Indonesia

PT. Interlab Sentra Solutions Indonesia

Tel (6221) 77840996

E-mail: office@issi-interlab.com

United Arab Emirates

Pasteur Central Labs Ltd

Telephone: +971-2-4467034

E-mail: pasteur@emirates.net.ae

www.pcl-uae.com